

**Variety developed to a fine point.**

The wide range of options offered by the SONICflex, will soon make it one of your most-used handpieces.

Choose the right tip for your needs, from the wide range of options.

**Numerous possibilities – main features and advantages for you:**

- More tips and more uses than any other airscaler
- Over 50 different tips
- For every indication
- Ensuring successful, minimally-invasive treatment




**SONICflex scaler tip no. 5**

Treatment: for removal of large areas of calculus  
Form: Universal



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS



0.571.5171  
1.005.8949

**SONICflex scaler tip no. 8**

Therapy: for removal of supragingival and subgingival calculus  
Form: Perio, extra-long



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS



0.571.5371  
1.006.1953

**SONICflex scaler tip no. 6**

Therapy: for removal of large areas of supragingival calculus  
Form: Sickle



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS




0.571.5181  
1.005.8950

**SONICflex scaler tip no. 7**

Therapy: for removal of large areas of supragingival calculus  
Form: Perio



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS



0.571.5191  
1.005.8951

**SONICflex scaler tip no. 48**

Therapy: removal of plaque, cleansing of fissures, cleansing of cavities and crown preparations  
Form: various attachments, single-use brushes



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS  
clean brush no. 1  
clean brush no. 2  
clean brush no. 3  
clean brush no. 4  
clean brush no. 5  
clean brush no. 6  
tip no. 48  
tip no. 48 A




0.571.0001  
1.006.2026  
0.571.0402  
0.571.0412  
0.571.0422  
0.571.0432  
0.571.0472  
0.571.0482  
0.571.0401  
1.006.1982

Prophylaxis

**SONICflex paro tip no. 60**

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; for the anterior region, smooth surfaces

Form: long straight tip - rounded surface



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7402**  
**1.006.1934**

**SONICflex paro tip no. 61**

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; from first/third quadrant to central incisor accessed buccally, then cleansing of front teeth, premolars and molars of the second/fourth quadrant, accessed palatally / lingually

Form: angled tip, **left** - rounded surface



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7412**  
**1.006.1935**

**SONICflex paro tip no. 62**

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; from second/fourth quadrant to central incisor accessed buccally, then cleansing of front teeth, premolars and molars of the first/third quadrant, accessed palatally / lingually"

Form: angled tip, **right** - rounded surface




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7422**  
**1.006.1936**

**SONICflex implant tip**

Therapy: effective removal of subgingival plaque and calculus around implants

Form: conical tips, single-use plastic base (PEEK)



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS  
single-use plastic base ( 30 pcs.)

**1.003.8167**  
**1.006.2027**  
**1.003.8168**

Periodontology

**SONICflex rootplaner tip no. 24**

Therapy: for difficult-to-reach furcations, enables access to approximal furcations, even in cases of unbroken lines of dentition

Form: contra-angled button, small, **left**, diamond coating D25



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.5621**  
**1.006.1957**

**SONICflex rootplaner tip no. 27**

Therapy: effective cleansing of smooth surfaces and concave areas

Form: extra-long perio, button, large; diamond coating D25



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.5651**  
**1.006.1963**

**SONICflex rootplaner tip no. 25**

Therapy: for difficult-to-reach furcations, enables access to approximal furcations, even in cases of unbroken lines of dentition

Form: contra-angled button, small, **right**; diamond coating D25



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.5631**  
**1.006.1959**

**SONICflex rootplaner tip no. 26**

Therapy: effective cleansing of smooth surfaces and concave areas

Form: universal, button, small; diamond coating D25



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.5641**  
**1.006.1961**

Periodontology



**SONICflex micro tip no. 30**

Therapy: minimally-invasive treatment, for incipient caries  
 Form: small hemisphere, **mesial**, diamond coating D46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.6811**  
**1.006.1969**

**SONICflex micro tip no. 33**

Therapy: minimally-invasive treatment, for incipient caries  
 Form: large hemisphere, **distal**, diamond coating D46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.6821**  
**1.006.1975**

**SONICflex micro tip no. 31**

Therapy: minimally-invasive treatment, for incipient caries  
 Form: small hemisphere, **distal**, diamond coating D46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.6801**  
**1.006.1971**

**SONICflex bevel tip no. 58**

Therapy: bevelling of approximal cavity edges  
 Form: **mesial**, diamond coating D46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7392**  
**1.006.1988**

**SONICflex micro tip no. 32**

Therapy: minimally-invasive treatment, for incipient caries  
 Form: large hemisphere, **mesial**, diamond coating D46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.6831**  
**1.006.1973**

**SONICflex bevel tip no. 59**

Therapy: bevelling of approximal cavity edges  
 Form: **distal**, diamond coating D46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS


**0.571.7382**  
**1.006.1990**

Minimally-invasive surgery



**SONICflex seal tip no. 45**

Treatment: minimally-invasive preparation of fissures, no-pain when opening of fissures  
 Form: conical, straight tip Ø 0.5 mm, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.000.8323  
 1.007.1506

**SONICflex cariex TC tip no. 71**

Therapy: controlled, layer-by-layer removal of carious tissues  
 Form: ball, Ø 1.0 mm, hard metal




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.000.7362  
 1.006.2002

**SONICflex cariex D tip no. 42**

Therapy: ultra-conservative, sonoabrasive caries excavation, fine  
 Form: ball, Ø 0.8 mm, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.000.7105  
 1.006.1980

**SONICflex cariex TC tip no. 72**

Therapy: controlled, layer-by-layer removal of carious tissues  
 Form: ball, Ø 1.4 mm, hard metal





Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.000.7363  
 1.006.2004

**SONICflex cariex D tip no. 43**

Therapy: ultra-conservative, sonoabrasive caries excavation, fine  
 Form: ball, Ø 1.2 mm, diamond coating D 64



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.000.7167  
 1.006.1981

Minimally-invasive surgery






**SONICflex prep gold tip no. 49, mesial**

Therapy: precise preparation of various bevel-angles around margins, after "rough preparation" with rotary instrument; this enables improved dental restoration

Form: **mesial**, tips with defined bevel-angles: lateral - 60° cervical - 75°, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7212**  
**1.006.1983**

**SONICflex prep ceram tip no. 52**

Treatment: precise preparation of various bevel-angles around margins, with gradual angle transitions from cervical to lateral areas

Form: **distal**, tips with defined bevel-edges, lateral - 60°, cervical - 75°, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7272**  
**1.006.1986**

**SONICflex prep gold tip no. 50 distal**

Therapy: precise preparation of various bevel-angles around margins, after "rough preparation" with rotary instruments; this enables improved dental restoration

Form: **distal**, tips with defined bevel-edges: lateral - 60°, cervical - 75°, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**0.571.7222**  
**1.006.1984**

**SONICflex prep ceram tip no. 51**

Treatment: precise preparation of various bevel-angles around margins, with gradual angle transitions from cervical to lateral areas

Form: **mesial**, tips with defined bevel-angles: lateral - 60°, cervical - 75°, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS


**0.571.7262**  
**1.006.1985**

Prosthetics



**SONICflex prep CAD / CAM tip no. 34**

Therapy: perfect adaptation of cavity preparations, to the specific requirements of CAD-CAM systems  
 Form: **mesial**, lateral angle 60°, cervical angle 90°, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**1.002.1984**  
**1.006.1977**

**SONICflex prep CAD / CAM tip no. 35**

Therapy: perfect adaptation of cavity preparations, to the specific requirements of CAD-CAM systems  
 Form: **distal**, lateral angle 60°, cervical angle 90°, diamond coating D 46





Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

**1.002.1986**  
**1.006.1979**

**SONICflex cem tip no. 12**

Therapy: gentle and secure placement of inlays  
 Form: rounded plastic cover



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS  
 CEM attachment



**0.571.5431**  
**1.006.2021**  
**0.571.7142**

**SONICflex crown prep tip no.97**

\*available from 16 May

Therapy: crown preparation - the use of rotary, crown preparation instruments, in combination with sonically-vibrated tips of similar form, is the perfect protocol for successful crown preparations.

Form: radius, diamond coating D46



Order number: compatible with SONICflex 2000 N / L and 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS



**1.008.6383**  
**1.008.6384**

**SONICflex crown prep tip no.98**

\*available as of 16 May

Therapy: crown preparation - the use of rotary, crown preparation instruments, in combination with sonically-vibrated tips of similar form, is the perfect protocol for successful crown preparations.

Form: radius, diamond coating D46



Order number: compatible with SONICflex 2000 N / L and 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S

**1.008.6385**  
**1.008.6386**

Prosthetics






**SONICflex retro tip no. 56**

Therapy: searching for/locating of root canal; significant work simplification during retrograde cavity preparation

Form: extremely narrow cylindrical tip, **left**, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

0.571.7322  
1.006.2040

**SONICflex retro tip no. 17**

Therapy: preparation of root canal, significantly simplified access, even with highly-convoluted roots

Form: **right**, small cylinder, diamond coating D 64




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

0.571.5581  
1.006.2035

**SONICflex retro tip no. 57**

Therapy: searching for/locating of root canal; significant work simplification during retrograde cavity preparation

Form: narrow cylindrical tip, **right**, diamond coating D 46




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

0.571.7332  
1.006.2041

**SONICflex retro tip no. 20 and no. 21**

Therapy: preparation of the undercut; correct, internal circumferential root canal preparation, with minimal bone defects

Form: pear-shaped, **small, left** (no. 20), **right** (no. 21), diamond coating D 64



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
matching for SONICflex 2008 / L, 2008 S / LS

0.571.5521  
1.006.2037




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
matching for SONICflex 2008 / L, 2008 S / LS

0.571.5561  
1.006.2038

**SONICflex retro tip no. 16**

Therapy: preparation of root canal, significantly simplified access, even with highly-convoluted roots

Form: **left**, small cylinder, diamond coated D 64



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

0.571.5541  
1.006.2036

**SONICflex retro tip no. 55**

Therapy: for roots, with apices highly-inclined towards buccal cavity

Form: diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

0.571.7342  
1.006.2039

Root tip resection



**SONICflex endo tip no. 66**

Therapy: shaping of access cavity, removal of protuberances  
 Form: button, **large**, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS


**1.000.5825**  
**1.006.1992**

**SONICflex endo tip no. 67**

Therapy: expansion of larger channels  
 Form: conical **125°**, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS



**1.000.5822**  
**1.006.1994**

**SONICflex endo tip no. 68**

Therapy: expansion of larger channels  
 Form: conical **112°**, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS


**1.000.5823**  
**1.006.1996**

**SONICflex endo tip no. 69**

Therapy: removal of residual root canal filling materials  
 Form: button, **small**, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS




**1.000.5827**  
**1.006.1998**

**SONICflex endo tip no. 70**

Therapy: for particularly long and wide canals  
 Form: conical **117°**, diamond coating D 25  
 anterior area is not diamond coated



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS



**1.000.5821**  
**1.006.2000**

**SONICflex tip Endo clean no. 96**

\*available from 16 May


Therapy: for flushing of root canal, for removal of debris and smear layer, as well as enhancing the performance of antibacterial rinsing solution  
 Form: 105° resp. 75° angle



Order number: compatible with SONICflex 2008 / L, 2008 S / LS

SONICflex Endo clean set A **1.007.1142**  
 SONICflex Endo clean tip A **1.008.5164**  
 compatible with SONICflex 2000 N / L and 2003 / L  
 SONICflex Endo clean tip **1.008.5163**  
 SONICflex Endo clean set **1.006.2045**

Endo clean needles for SONICflex 2000 N / L, 2003 and 2008  
 SONICflex Endo clean 015 white **1.006.2042**  
 SONICflex Endo clean 020 yellow **1.006.2043**  
 SONICflex Endo clean 025 red **1.006.2044**




Endodontics



**SONICflex bone tip no. 80**

Therapy: sectioning of Sharpey's fibres; extension of alveolae; removal of root residues; gentle extraction before implantation; bone splitting

Form: square, cutting




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**1.004.3875**  
**1.006.2006**

**SONICflex bone tip no. 84**

Therapy: removal of bone blocks, osteotomies (distraction, bone splitting); if applicable, dealing with extremely thick, substantia corticalis

Form: saw axial




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**1.004.3879**  
**1.006.2014**

**SONICflex bone tip no. 81**

Therapy: preparation of bone window for sinus lift; with very massive bones, if applicable, rotating and sectioning of substantia corticalis, plus removal of remaining bone using the diamond-coated tip

Form: ball, large, diamond coating D 46



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**1.004.3876**  
**1.006.2008**

**SONICflex bone tip no. 85**

Therapy: the disc-shaped form, makes it ideal for safe repositioning of the Schneider membrane, during a sinus lift

Form: elephant foot



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**1.006.0645**  
**1.007.1624**

**SONICflex bone tip no. 83**

Therapy: removal of bone blocks, osteotomies (distraction, bone splitting); if applicable, dealing with extremely thick, substantia corticalis

Form: saw sagittal




Order number: compatible with SONICflex 2000 N / L, 2003 / L  
compatible with SONICflex 2008 / L, 2008 S / LS

**1.004.3878**  
**1.006.2012**

Surgical

**SONICflex bone tip no. 86**

Therapy: ideal for easy and quick harvesting, of greater quantities of bone, especially useful for bone augmentation  
 Form: bone harvester



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS

1.006.0741  
 1.007.1625

**SONICflex bone tip no. 87**

Therapy: ideal for precise, extremely narrow incisions, when cutting ligaments  
 Form: saw axial, replaceable saw blades  
 Saw blade thickness: 0.15 mm



Order number: compatible with SONICflex 2000 N / L, 2003 / L  
 compatible with SONICflex 2008 / L, 2008 S / LS  
 Saw blade ( 5 pcs.)

1.006.1406  
 1.007.1627  
 1.006.1405

**Accessories:**



Order number: torque wrench

1.000.4887



Order number: steri cassette

0.411.9101



Order number: key

0.411.0892



Order number: adapter for supply of liquid medicaments, for SONICflex 2003  
 adapter for supply of liquid medicaments for SONICflex 2008

1.000.6296  
 1.006.5966

Surgical

